



Product data sheet REF 68.83.44A


25.08.2026

General data			
article number	68.83.44A	Mikro screw self-drill green 1.5x3.5mm P/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	46638	GTIN	04046826541699
Risk category	N.A. (EU MDD)	UMDNS-Code	
Classification	NA (EU MDD)	Base UDI-DI	4046826;B001;01;P8
Net weight (kg)	0,000	Operating instructions	G687004USA
Latex free	yes	Preparation instructions	MPK06



Medicon media

Complete picture	Detail images
------------------	---------------



size	1.5x3.5mm	Coating	
------	-----------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))