






## Product data sheet REF 68.84.07A

11.12.2025

General data			
article number	68.84.07A	Ti Mini locking plate 6 hole 60° twisted lf 1.0mm fracture	
Unit	piece	Customs code (HS code)	90189084
GMDN	46642	GTIN	04046826434175
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
size	60° 1.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		