






## Product data sheet REF 68.84.08

11.12.2025

| General data  |                                   |  |                |
|---|-----------------------------------|--|----------------|
| article number  | 68.84.08                          | Ti Mini locking plate 6 hole 60° twisted rt 1.0mm fracture |                |
| Unit  | piece                             | Customs code (HS code)                                     | 90189084       |
| GMDN  | 46642                             | GTIN   | 04046826434182 |
| Risk category   | 2b (EU MDD)                       | UMDNS-Code   | 13-050         |
| Classification  | 8 (EU MDD)                        | Base UDI-DI  |                |
| Net weight (kg)   | 0,000                             | Operating instructions                                     | G689480        |
| Latex free  | yes                               | Preparation instructions                                   | MPK06          |
|    |                                   |  |                |
| Medicon media   |                                   |  |                |
| Complete picture  |                                   | Detail images  |                |
|   |                                   |  |                |
| size  | 60° 1.0mm                         | Coating  |                |
| Storage conditions  |                                   |  |                |
| No special requirements   |                                   |  |                |
| Materials   |                                   |  |                |
|   | Material number (DIN designation) |  |                |
| Main part   | 3.7035 (ASTM_Grade_2)             |  |                |