







Product data sheet REF 68.84.21A

11.12.2025

| General data | | | |
|---|-----------------------------------|--|----------------|
| article number | 68.84.21A | Ti Mini locking plate 5x10 hole ang lf 2.0mm fract | |
| Unit | piece | Customs code (HS code) | 90211090 |
| GMDN | 46642 | GTIN | 04046826345310 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 13-050 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,000 | Operating instructions | G689480 |
| Latex free | yes | Preparation instructions | MPK06 |
|    | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 2.0mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.7035 (ASTM_Grade_2) | | |