







Product data sheet REF 68.84.22A

11.12.2025

General data			
article number	68.84.22A	Ti Mini locking plate 6 holes ang 2.0mm fracture	
Unit	piece	Customs code (HS code)	90211090
GMDN	46642	GTIN	04046826345327
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	2.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		