







## Product data sheet REF 68.84.42A

11.12.2025

General data			
article number	68.84.42A	Ti Mini locking plate 6x20 hole ang lf 2.0mm fract	
Unit	piece	Customs code (HS code)	90211090
GMDN	46642	GTIN	04046826427634
Risk category	2b (EU MDD)	UMDNS-Code	13-050
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
			
size	2.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7035 (ASTM_Grade_2)		