



Product data sheet REF 68.85.43A

11.12.2025

| General data | | | |
|---|---|---|----------------|
| article number | 68.85.43A | Ti Micro cross slot screw 1.5x3.5mm grey self-tapp p/10 | |
| Unit | box/10 | Customs code (HS code) | 90211090 |
| GMDN | 61837 | GTIN | 04046826345914 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 16-101 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,001 | Operating instructions | G687004 |
| Latex free | yes | Preparation instructions | MPK06 |
|    | | | |
| Medicon media | 451.08.13 (brochure intraoral distraction devices f upper a lower jaws) 453.17.00_DEU (Le Fort I Plate acc to Rustemeyer german INTRO) 453.17.00_ENG (Le Fort I Plate acc to Rustemeyer english INTRO) 453.17.00_ESP (Le Fort I Plate acc to Rustemeyer spanish INTRO) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 1.5x3.5mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI)) | | |