



## Product data sheet REF 68.85.57A

11.12.2025

| General data    |             |                                                        |                |
|-----------------|-------------|--------------------------------------------------------|----------------|
| article number  | 68.85.57A   | Ti Micro cross slot screw 1.5x17mm grey self-tapp p/10 |                |
| Unit            | box/10      | Customs code (HS code)                                 | 90211090       |
| GMDN            | 61837       | GTIN                                                   | 04046826404017 |
| Risk category   | 2b (EU MDD) | UMDNS-Code                                             | 16-101         |
| Classification  | 8 (EU MDD)  | Base UDI-DI                                            |                |
| Net weight (kg) | 0,000       | Operating instructions                                 | G687004        |
| Latex free      | yes         | Preparation instructions                               | MPK06          |



Medicon media

## Complete picture

## Detail images



|      |          |         |  |
|------|----------|---------|--|
| size | 1.5x17mm | Coating |  |
|------|----------|---------|--|

## Storage conditions

No special requirements

## Materials

|           | Material number (DIN designation)   |
|-----------|-------------------------------------|
| Main part | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI)) |