






## Product data sheet REF 68.86.44A

11.12.2025

| General data  |                                   |                                 |                |
|---|-----------------------------------|---------------------------------|----------------|
| article number  | 68.86.44A                         | Ti Micro Le Fort I plate lf 5mm |                |
| Unit  | piece                             | Customs code (HS code)          | 90211090       |
| GMDN  | 46642                             | GTIN                            | 04046826427740 |
| Risk category   | 2b (EU MDD)                       | UMDNS-Code                      | 13-050         |
| Classification  | 8 (EU MDD)                        | Base UDI-DI                     |                |
| Net weight (kg)   | 0,000                             | Operating instructions          | G687004        |
| Latex free  | yes                               | Preparation instructions        | MPK06          |
|    |                                   |                                 |                |
| Medicon media   |                                   |                                 |                |
| Complete picture  |                                   | Detail images                   |                |
|   |                                   |                                 |                |
| size  | 5mm                               | Coating                         |                |
| Storage conditions  |                                   |                                 |                |
| No special requirements   |                                   |                                 |                |
| Materials   |                                   |                                 |                |
|   | Material number (DIN designation) |                                 |                |
| Main part   | 3.7035 (ASTM_Grade_2)             |                                 |                |