







## Product data sheet REF 68.87.78A

11.12.2025

| General data  |                                   |   |                |
|---|-----------------------------------|---|----------------|
| article number  | 68.87.78A                         | Ti Micro double Y plate short Rigid p/5 |                |
| Unit  | box/5                             | Customs code (HS code)                  | 90211090       |
| GMDN  | 46642                             | GTIN                                    | 04046826417437 |
| Risk category   | 2b (EU MDD)                       | UMDNS-Code                              | 13-050         |
| Classification  | 8 (EU MDD)                        | Base UDI-DI                             |                |
| Net weight (kg)   | 0,000                             | Operating instructions                  | G687004        |
| Latex free  | yes                               | Preparation instructions                | MPK06          |
|    |                                   |   |                |
| Medicon media   |                                   |   |                |
| Complete picture  |                                   | Detail images                           |                |
|   |                                   |   |                |
| Storage conditions  |                                   |   |                |
| No special requirements   |                                   |   |                |
| Materials   |                                   |   |                |
|   | Material number (DIN designation) |   |                |
| Main part   | 3.7035 (ASTM_Grade_2)             |   |                |