



Product data sheet REF 68.88.03

11.12.2025

| General data | | | |
|---|---|--|----------------|
| article number | 68.88.03 | mandib distract dev horiz 20mm f lf side | |
| Unit | piece | Customs code (HS code) | 90211090 |
| GMDN | 47635 | GTIN | 04046826126032 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 11-291 |
| Classification | 8 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,042 | Operating instructions | G688802 |
| Latex free | yes | Preparation instructions | MPK06 |
|    | | | |
| Medicon media | 451.08.13 (brochure intraoral distraction devices f upper a lower jaws) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | 20mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 3.7165 (ASTM_Grade_5_TiAl6V4_(ELI)) | | |
| Screw | 1.4441 (X2CrNiMo18-15-3) | | |
| Others | 3.7035 (ASTM_Grade_2) | | |