



Product data sheet REF 68.88.63

11.12.2025

General data			
article number	68.88.63	transantral maxilla distractor 25mm f lf side	
Unit	piece	Customs code (HS code)	90211090
GMDN	47635	GTIN	04046826304171
Risk category	2b (EU MDD)	UMDNS-Code	11-291
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G688862
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	451.08.13 (brochure intraoral distraction devices f upper a lower jaws) 451.08.22 (brochure transantral maxilla distraction devices)		
Complete picture		Detail images	
 			
size	25mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		
Screw	1.4441 (X2CrNiMo18-15-3)		
Others	3.7035 (ASTM_Grade_2)		