







## Product data sheet REF 68.90.24

11.12.2025

| General data  |                                   |  |                |
|---|-----------------------------------|--|----------------|
| article number  | 68.90.24                          | Ti Mini compression plate 4 holes long |                |
| Unit  | piece                             | Customs code (HS code)                 | 90211090       |
| GMDN  | 46642                             | GTIN                                   | 04046826126247 |
| Risk category   | 2b (EU MDD)                       | UMDNS-Code                             | 13-050         |
| Classification  | 8 (EU MDD)                        | Base UDI-DI                            |                |
| Net weight (kg)   | 0,010                             | Operating instructions                 | G687004        |
| Latex free  | yes                               | Preparation instructions               | MPK06          |
|    |                                   |  |                |
| Medicon media   |                                   |  |                |
| Complete picture  |                                   | Detail images                          |                |
|   |                                   |  |                |
| Storage conditions  |                                   |  |                |
| No special requirements   |                                   |  |                |
| Materials   |                                   |  |                |
|   | Material number (DIN designation) |  |                |
| Main part   | 3.7035 (ASTM_Grade_2)             |  |                |