






Product data sheet REF 68.92.18A

11.12.2025

General data			
article number	68.92.18A	Ti CMS Emerg cross slot screw 1.4x3mm red self-tapping p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826378615
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689000
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	453.12.00 (Epiplating)		
Complete picture		Detail images	
size	1.4x3mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		