






## Product data sheet REF 68.92.19A

11.12.2025

General data			
article number	68.92.19A	Ti CMS Emerg cross slot screw 1.4x5mm red self-tapping p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826378639
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689000
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	453.12.00 (Epiplating )		
Complete picture		Detail images	
size	1.4x5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		