



Product data sheet REF 68.94.00A

11.12.2025

General data			
article number	68.94.00A	Ti Mini cross slot screw 2.0x4.0mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826349448
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
			
Medicon media	451.09.02 (brochure bimaxillary surgery) 453.12.00 (Epiplating)		
Complete picture		Detail images	
size	2.0x4.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		