



Product data sheet REF 68.94.03A

11.12.2025

General data			
article number	68.94.03A	Ti Mini cross slot screw 2.0x9.0mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826349509
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	451.08.13 (brochure intraoral distraction devices f upper a lower jaws) 451.09.02 (brochure bimaxillary surgery) 453.06.00_DEU (Hybrid Plate System german Brochure) 453.06.00_ENG (Hybrid Plate System - Reconstruction of the mandible with a microvascular fibular bone graft ENG Brochure) 453.06.00_ESP (Hybrid Plate System - Reconstruction of the mandible with a microvascular fibular bone graft ESP Brochure)		
Complete picture		Detail images	
size	2.0x9.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		