



## Product data sheet REF 68.94.05A

11.12.2025

General data			
article number	68.94.05A	Ti Mini cross slot screw 2.0x13mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826349547
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G687004
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media	451.09.02 (brochure bimaxillary surgery )		
Complete picture		Detail images	
size	2.0x13mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		