



Product data sheet REF 68.94.10L

11.12.2025

General data			
article number	68.94.10L	Ti Mini Emerg cross slot screw 2.3x7.0mm red self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826470845
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,006	Operating instructions	G687004L
Latex free	yes	Preparation instructions	MPK06
			
Medicon media			
Complete picture		Detail images	
size	2.3x7.0mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		