



Product data sheet REF 68.94.53A

11.12.2025

General data			
article number	68.94.53A	Ti cross slot screw 2.4x13mm grey self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826349974
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
size	2.4x13mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		