



Product data sheet REF 68.94.90A

11.12.2025

General data			
article number	68.94.90A	Ti Mini cross slot screw lock 2.0x7mm gold self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826350420
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,001	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
			
Medicon media			
Complete picture		Detail images	
			
size	2.0x7mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		