



Product data sheet REF 68.94.93A

11.12.2025

General data			
article number	68.94.93A	Ti Mini cross slot screw lock 2.0x15mm gold self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826350482
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06



Medicon media

Complete picture

Detail images



size	2.0x15mm	Coating	
------	----------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))