






## Product data sheet REF 68.95.34

11.12.2025

General data			
article number	68.95.34	plate bending lever 15cm f syst 2.4	
Unit	piece	Customs code (HS code)	90189084
GMDN	44795	GTIN	04046826127503
Risk category	1 (EU MDR)	UMDNS-Code	11-508
Classification	1 (EU MDR)	Base UDI-DI	4046826;N001;06;WT
Net weight (kg)	0,080	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
			
size	15cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		