



Product data sheet REF 68.95.72

11.12.2025

General data			
article number	68.95.72	twist drill Ø 1.8x25mm 50mm Stryker	
Unit	piece	Customs code (HS code)	90184910
GMDN	10520	GTIN	04046826336806
Risk category	2a (EU MDD)	UMDNS-Code	11-341
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media			
Complete picture		Detail images	
			
size	Ø 1.8x25mm 50mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		