




Product data sheet REF 68.98.25A

11.12.2025

General data			
article number	68.98.25A	Ti Star Grip screw locking 2.4x15mm gold self-tapp p/10	
Unit	box/10	Customs code (HS code)	90211090
GMDN	61837	GTIN	04046826351595
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,005	Operating instructions	G689480
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
		  Selbstschneidend Self-tapping Autofrescante Autofrauseuse Autoflettante	
size	2.4x15mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		