



## Product data sheet REF 68.99.66

11.12.2025

General data			
article number	68.99.66	Maxpander anchoring screw self-dr f 2.0 P/4	
Unit	box/4	Customs code (HS code)	90211090
GMDN	64609	GTIN	04046826466886
Risk category	2b (EU MDD)	UMDNS-Code	16-101
Classification	8 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,008	Operating instructions	G686690
Latex free	yes	Preparation instructions	MPK06
  			
Medicon media			
Complete picture		Detail images	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.7165 (ASTM_Grade_5_TiAl6V4_(ELI))		