





Product data sheet REF 681.00.34

11.12.2025

General data			
article number	681.00.34	Mixer clamp 90° Ø 5mm double action insert	
Unit	piece	Customs code (HS code)	90189084
GMDN	33596	GTIN	04046826306724
Risk category	2b (EU MDD)	UMDNS-Code	11-502
Classification	9 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,026	Operating instructions	G6810101
Latex free	yes	Preparation instructions	MPK04
			
Medicon media			
Complete picture		Detail images	
			
size	90° Ø 5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Others	1.4310 (X10CrNi18-8)		