





## Product data sheet REF 681.01.04

11.12.2025

| General data   |  |   |  |
|--|--|---|--|
| article number   | 681.01.04                                    | hook scrs Ø 5mm single action w/o ratchet |  |
| Unit   | piece  | Customs code (HS code)                    | 90189084   |
| GMDN   | 33596  | GTIN                                      | 04046826306083   |
| Risk category  | 2b (EU MDD)                                  | UMDNS-Code                                | 13-484   |
| Classification   | 9 (EU MDD)                                   | Base UDI-DI                               |  |
| Net weight (kg)  | 0,126  | Operating instructions                    | G6810101   |
| Latex free   | yes  | Preparation instructions                  | MPK04  |
|   |  |   |  |
| Medicon media  |  |   |  |
| Complete picture   |  | Detail images                             |  |
|  |  |   |  |
| size   | Ø 5mm  | Coating                                   | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions   |  |   |  |
| No special requirements  |  |   |  |
| Materials  |  |   |  |
|  | Material number (DIN designation)            |   |  |
| <b>Main part</b>   | 1.4021 (X20Cr13)                             |   |  |
| <b>Screw</b>   | 1.4305 (X8CrNiS18-9)                         |   |  |
| <b>Handle</b>  | PA6.6GF30 (Polyamid mit 30% Glasfaseranteil) |   |  |
| <b>Others</b>  | PPSU (Polyphenylsulfon)                      |   |  |