







Product data sheet REF 681.03.30

11.12.2025

| General data | | | |
|--|--|--|--|
| article number | 681.03.30 | dolphin grasp fcps Ø 5mm double action w ratchet | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 33596 | GTIN | 04046826306557 |
| Risk category | 2b (EU MDD) | UMDNS-Code | 11-502 |
| Classification | 9 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,221 | Operating instructions | G6810101 |
| Latex free | yes | Preparation instructions | MPK04 |
|  | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
|    | | | |
| size | Ø 5mm | Coating | Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Screw | 1.4305 (X8CrNiS18-9) | | |
| Handle | PA6.6GF30 (Polyamid mit 30% Glasfaseranteil) | | |
| Others | PPSU (Polyphenylsulfon) | | |