







Product data sheet REF 681.03.34

11.12.2025

General data			
article number	681.03.34	Mixer clamp 90° Ø 5mm double action w ratchet	
Unit	piece	Customs code (HS code)	90189084
GMDN	33596	GTIN	04046826306717
Risk category	2b (EU MDD)	UMDNS-Code	11-502
Classification	9 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,129	Operating instructions	G6810101
Latex free	yes	Preparation instructions	MPK04
			
Medicon media			
Complete picture		Detail images	
  			
size	90° Ø 5mm	Coating	Artificially created galvanic gold layer (max 0.04mm) on parts of the surface of instruments
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		
Screw	1.4305 (X8CrNiS18-9)		
Handle	PA6.6GF30 (Polyamid mit 30% Glasfaseranteil)		
Others	PPSU (Polyphenylsulfon)		