



Product data sheet REF 69.02.10

11.12.2025

General data			
article number	69.02.10	implant tray Mini Fracture Locking Syst 2.0	
Unit	piece	Customs code (HS code)	90189084
GMDN	12143	GTIN	04046826338725
Risk category	1 (EU MDR)	UMDNS-Code	12-142
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;03;X8
Net weight (kg)	16,400	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
			
size		Coating	Anodizing is the process of enriching titanium with a hard anodized layer on aluminum components.
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.3547 (AlMg4,5Mn)		
Screw	1.4305 (X8CrNiS18-9)		
Others	PPSU (Polyphenylsulfon)		