



## Product data sheet REF 69.02.80

11.12.2025

General data			
article number	69.02.80	instr tray Micro/Mini Syst 1.5/2.0	
Unit	piece	Customs code (HS code)	90189084
GMDN	12143	GTIN	04046826371944
Risk category	1 (EU MDR)	UMDNS-Code	12-142
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;03;X8
Net weight (kg)	2,054	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
 			
size		Coating	Anodizing is the process of enriching titanium with a hard anodized layer on aluminum components.
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.3547 (AlMg4,5Mn)		
Screw	1.4305 (X8CrNiS18-9)		
Others	PPSU (Polyphenylsulfon)		