






Product data sheet REF 69.03.01

11.12.2025

General data			
article number	69.03.01	outer tray w 2 storage shelves f impl trays 1.2 to 2.4 fract	
Unit	piece	Customs code (HS code)	90189084
GMDN	12143	GTIN	04046826353391
Risk category	1 (EU MDR)	UMDNS-Code	
Classification	1 (EU MDR)	Base UDI-DI	4046826;N003;03;X8
Net weight (kg)	1,050	Operating instructions	GATDN00X
Latex free	yes	Preparation instructions	
 			
Medicon media			
Complete picture		Detail images	
			
size		Coating	Anodizing is the process of enriching titanium with a hard anodized layer on aluminum components.
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	3.3547 (AlMg4,5Mn)		
Screw	1.4305 (X8CrNiS18-9)		