





Product data sheet REF 70.35.03

11.12.2025

| General data | | | |
|---|-----------------------------------|-----------------------------------|-------------------------|
| article number | 70.35.03 | Olivecrona dissector 19.5cm 4+5mm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 39610 | GTIN | 04046826493042 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 15-249 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;04;aj;5X |
| Net weight (kg) | 0,020 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | | | |
| Complete picture | | Detail images | |
| | | | |
| size | 19.5cm 4+5mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |