



Product data sheet REF 71.32.08

11.12.2025

| General data | | | |
|-----------------|-------------|------------------------------|-------------------------|
| article number | 71.32.08 | Vogel adenoid curette 21.5cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 33449 | GTIN | 04046826129682 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 16-583 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;02;ag;37 |
| Net weight (kg) | 0,047 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK12-1 |



| | |
|---------------|----------------------------|
| Medicon media | 450.00.37 (catalogue ENT) |
|---------------|----------------------------|

Complete picture

Detail images



| | | | |
|------|--------|---------|--|
| size | 21.5cm | Coating | |
|------|--------|---------|--|

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|---------------------|-----------------------------------|
| Main part Handle | 1.4021 (X20Cr13) |
| | 1.4301 (X5CrNi18-10) |