



Product data sheet REF 71.36.18

11.12.2025

| General data | | | |
|--|-----------------------------------|-------------------------------------|-------------------------|
| article number | 71.36.18 | Barnhill adenoid curette 22cm fig 4 | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 33449 | GTIN | 04046826129873 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 16-583 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R001;02;ag;37 |
| Net weight (kg) | 0,060 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | MPK12-1 |
|   | | | |
| Medicon media | 450.00.37 (catalogue ENT) | | |
| Complete picture | | Detail images | |
|   | | | |
| size | 22cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Handle | 1.4301 (X5CrNi18-10) | | |