



## Product data sheet REF 710.30.13

11.12.2025

General data			
article number	710.30.13	drill Ø 2.5mm 110mm f quick coupling	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826357740
Risk category	2a (EU MDD)	UMDNS-Code	15-204
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	Ø 2.5mm 110mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4112 (X90CrMoV18)		
Others	1.4021 (X20Cr13)		