





## Product data sheet REF 710.30.21

11.12.2025

General data			
article number	710.30.21	tap 110mm f 3.5mm f quick coupling	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826357825
Risk category	2a (EU MDD)	UMDNS-Code	11-331
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	110mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4034 (X46Cr13)		