



## Product data sheet REF 711.30.20

11.12.2025

General data			
article number	711.30.20	drill Ø 4.5mm 145mm f quick coupling	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826361471
Risk category	2a (EU MDD)	UMDNS-Code	15-204
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	Ø 4.5mm 145mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		