



Product data sheet REF 711.30.29

11.12.2025

General data			
article number	711.30.29	tap 125mm f 4.5mm f quick coupling	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826361563
Risk category	2a (EU MDD)	UMDNS-Code	11-331
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,000	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
			
size	125mm 4.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4034 (X46Cr13)		