



Product data sheet REF 711.30.31

11.12.2025

General data			
article number	711.30.31	tap 175mm f 4.5mm calibrated	
Unit	piece	Customs code (HS code)	90189084
GMDN	32390	GTIN	04046826361587
Risk category	2a (EU MDD)	UMDNS-Code	11-331
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,012	Operating instructions	G930100
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.63 (catalogue Orthopaedics)		
Complete picture		Detail images	
 			
size	175mm 4.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4034 (X46Cr13)		