



## Product data sheet REF 712.35.52

11.12.2025

General data			
article number	712.35.52	countersink 160mm f screws cannulated Ø 3.5mm	
Unit	piece	Customs code (HS code)	90189084
GMDN	16462	GTIN	04046826376222
Risk category	2a (EU MDD)	UMDNS-Code	11-331
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,017	Operating instructions	
Latex free	yes	Preparation instructions	MPK02
			
Medicon media	450.00.63 (catalogue Orthopaedics )		
Complete picture		Detail images	
			
size	160mm Ø 3.5mm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		