



Product data sheet REF 713.30.01

11.12.2025

| General data | | | |
|--|-------------------------------------|---|----------------|
| article number | 713.30.01 | drill three-fluted Ø 3.2mm 230mm f quick coupling | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 32390 | GTIN | 04046826370527 |
| Risk category | 2a (EU MDD) | UMDNS-Code | 15-204 |
| Classification | 6 (EU MDD) | Base UDI-DI | |
| Net weight (kg) | 0,016 | Operating instructions | G930100 |
| Latex free | yes | Preparation instructions | MPK02 |
|  | | | |
| Medicon media | 450.00.63 (catalogue Orthopaedics) | | |
| Complete picture | | Detail images | |
|  | | | |
| size | Ø 3.2mm 230mm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |