



Product data sheet REF 72.25.35

11.12.2025

General data			
article number	72.25.35	Henke tonsil dissector d/e 24.5cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	12255	GTIN	04046826131166
Risk category	1r (EU MDR)	UMDNS-Code	14-070
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;03;ai;4W
Net weight (kg)	0,048	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK12-1



Medicon media	450.00.37 (catalogue ENT)
---------------	----------------------------

Complete picture

Detail images



size	24.5cm	Coating	
------	--------	---------	--

Storage conditions

No special requirements

Materials

	Material number (DIN designation)
Main part	1.4021 (X20Cr13)
Handle	1.4301 (X5CrNi18-10)