







Product data sheet REF 72.25.38

11.12.2025

| General data | | | |
|---|-----------------------------------|--|-------------------------|
| article number | 72.25.38 | Henke tonsil dissector large d/e 24cm | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 12255 | GTIN | 04046826131180 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 14-070 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R004;04;aj;5X |
| Net weight (kg) | 0,046 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |
|   | | | |
| Medicon media | 450.00.37 (catalogue ENT) | | |
| Complete picture | | Detail images | |
| | |   | |
| | | | |
| size | 24cm | Coating | |
| Storage conditions | | | |
| No special requirements | | | |
| Materials | | | |
| | Material number (DIN designation) | | |
| Main part | 1.4021 (X20Cr13) | | |
| Handle | 1.4301 (X5CrNi18-10) | | |