







## Product data sheet REF 72.40.18

11.12.2025

General data			
article number	72.40.18	White tonsil seizing fcps 18cm	
Unit	piece	Customs code (HS code)	90189084
GMDN	62468	GTIN	04046826131241
Risk category	1r (EU MDR)	UMDNS-Code	14-070
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;02;ah;3V
Net weight (kg)	0,051	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK23-1
 			
Medicon media	450.00.37 (catalogue ENT )		
Complete picture		Detail images	
 			
size	18cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4021 (X20Cr13)		