






Product data sheet REF 72.50.95

11.12.2025

General data			
article number	72.50.95	wire loops soldered p/100	
Unit	box/100	Customs code (HS code)	90189084
GMDN	13632	GTIN	04046826131487
Risk category	1r (EU MDR)	UMDNS-Code	13-632
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R002;06;ah;59
Net weight (kg)	0,015	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	MPK01-1
 			
Medicon media	450.00.37 (catalogue ENT)		
Complete picture		Detail images	
			
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4310 (X10CrNi18-8)		