






Product data sheet REF 72.95.19

11.12.2025

General data			
article number	72.95.19	tonsil cannula Ø 0.8mm 12cm 90° ang LL	
Unit	piece	Customs code (HS code)	90183900
GMDN	60909	GTIN	04046826131647
Risk category	2a (EU MDD)	UMDNS-Code	15-206
Classification	6 (EU MDD)	Base UDI-DI	
Net weight (kg)	0,024	Operating instructions	
Latex free	yes	Preparation instructions	MPK03
			
Medicon media	450.00.37 (catalogue ENT)		
Complete picture		Detail images	
 			
size	Ø 0.8mm 12cm 90°	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part	1.4301 (X5CrNi18-10)		
Handle	2.0401 (CuZn39Pb3)		