








## Product data sheet REF 73.17.03

11.12.2025

General data			
article number	73.17.03	laryng cup fcps Ø 2mm 23cm rt cvd	
Unit	piece	Customs code (HS code)	90189084
GMDN	62467	GTIN	04046826132194
Risk category	1r (EU MDR)	UMDNS-Code	16-270
Classification	6 (EU MDR)	Base UDI-DI	4046826;B;R003;06;ai;5X
Net weight (kg)	0,040	Operating instructions	GATDR00X
Latex free	yes	Preparation instructions	
 			
Medicon media	450.00.37 (catalogue ENT )		
Complete picture		Detail images	
  			
size	Ø 2mm 23cm	Coating	
Storage conditions			
No special requirements			
Materials			
	Material number (DIN designation)		
Main part Screw	1.4021 (X20Cr13) 1.4305 (X8CrNiS18-9)		