



Product data sheet REF 73.17.04

11.12.2025

| General data | | | |
|-----------------|-------------|-----------------------------------|-------------------------|
| article number | 73.17.04 | laryng cup fcps Ø 2mm 23cm lf cvd | |
| Unit | piece | Customs code (HS code) | 90189084 |
| GMDN | 62467 | GTIN | 04046826132200 |
| Risk category | 1r (EU MDR) | UMDNS-Code | 16-270 |
| Classification | 6 (EU MDR) | Base UDI-DI | 4046826;B;R002;02;ah;3V |
| Net weight (kg) | 0,040 | Operating instructions | GATDR00X |
| Latex free | yes | Preparation instructions | |



Medicon media 450.00.37 (catalogue ENT)

Complete picture

Detail images



size Ø 2mm 23cm

Coating

Storage conditions

No special requirements

Materials

| | Material number (DIN designation) |
|--------------------|-----------------------------------|
| Main part Screw | 1.4021 (X20Cr13) |
| | 1.4305 (X8CrNiS18-9) |